

“ For Hygiene ”
BIORAD[®]

necm
kimya

Data Sheet

READY HP WIPE



Disinfectant wipes for surfaces and medical devices

- ✓ Disinfection for surfaces
- ✓ High level disinfectant wipes
- ✓ Disinfectant wipes based on hydrogen peroxide

PRODUCT SPECIFICATIONS

It is an anticorrosive product based on hydrogen peroxide, which has a pleasant smell and can be used for high-level disinfection of the surfaces of medical devices. Its formula is based on hydrogen peroxide and is a disinfectant that does not contain aldehyde, alcohol and phosphate.

FIELD OF USE

It is used to clean and disinfect all horizontal and vertical surfaces:

- surfaces made of: PVC, glass, PP, HDPE, acrylic glass, plexiglass, epoxy floors, walls, doors, patient beds.
- surfaces in health facilities, such as: medical offices, operating rooms, intensive therapy rooms, pediatrics, neonatology departments, dental offices, ENT offices, pneumophysiology departments, endoscopy offices.
- Homes for the elderly, schools, kindergartens, buses, ambulances, the pharmaceutical and cosmetics industry and other institutional spaces
- medical devices: ultrasound probes, stethoscopes, incubators, intravenous supports, microscopes, instruments that come into contact with the patient's skin.

APPLICATION METHOD

READY HP WIPE is packed in a dose of 100, 120, 300 napkins and is suitable for the disinfection of small surfaces and medical devices, by wiping. An action time of 5 minutes is necessary to achieve the effectiveness according to the desired spectrum.

After disinfection with the READY HP WIPE product, rinsing is not necessary.

Manufacturer: **Necm Kimya**

Address: Battalgazi Mah. Şark Cad. Aytop,

Gülöncaşı Sk. I Blok No:14/1 Sultanbeyli/İSTANBUL

Tel: +90 216 669 06 66

E-mail: info@necmkimya.com.tr

ACTIVITY / TESTING METHODS / ACTION TIME

Bactericidal		Minutes
E. hirae, S. Aureus, E. Coli (clean conditions)	EN 13697	1 min
E. hirae, S. Aureus, P. aeruginosa, E. Coli (dirty conditions)	EN 13697	5 min
E. hirae, P. aeruginosa, S. Aureus, C. albicans (clean conditions)	EN 16615	5 min
E. hirae, P. aeruginosa, S. Aureus (clean and dirty conditions)	EN 13727	1 min
Yeasticidal		
C. albicans, A. brasiliensis (clean and dirty conditions)	EN 13624	1 min
E. hirae, P. aeruginosa, S. Aureus, C. albicans (clean conditions)	EN 16615	5 min
Fungicide		
C. albicans, A. brasiliensis (clean and dirty conditions)	EN 13624	1 min
Mycobactericidal / Tuberculocidal		
M. terrae, M. avium (clean and dirty conditions)	EN 14348	1 min
Virucidal		
Adenovirus, Poliovirus, MNV, MVA (clean conditions)	EN 14476	5 min
Sporicidal		
B. subtilis (clean and dirty conditions)	EN 13704	1 min

MICROBIOLOGY

READY HP WIPE has bactericidal, fungicidal, mycobactericidal, tuberculocidal, virucidal, sporicidal activity

CHEMICAL-PHYSICAL DATA

➤ pH: 5,5 - 7

CHEMICAL COMPOSITION

Hydrogen peroxide 6 %- CAS 7722-84-1; CE 231-765-0

“ For Hygiene ”
BIORAD[®]



Data Sheet

READY HP WIPE

SHELF LIFE : 2 (two) years

PERIOD OF ACTIVITY AFTER OPENING THE COVER: 75 (seventy five) days

PACKAGING METHODS:

- Dose with 300 napkins 20x16 cm
- Dose with 120 napkins 14x16 cm
- Dose with 100 napkins 14x16 cm

CONDITIONS OF STORAGE AND DELIVERY OF PRODUCTS

- For the storage and preservation of products, the ambient temperature must be between 0 and 50 °C, and the humidity must be between 10% and 80% rH.
- It is sufficient to store the products above ground so that they are not exposed to direct sunlight and to prevent contact with water.
- A maximum of 7 (seven) boxes can be placed on top of each other both during storage and during loading, in order to preserve and not damage the packaging of the products.

GENERAL WARNINGS

- It is not classified as toxic, corrosive, mutagenic, irritating to the skin and does not cause respiratory tract sensitivity.
- In case of direct contact with the eyes, rinse immediately with plenty of water.
- It should not be used for cleaning hands, face and food.
- It should not be kept out of the reach of children.

SAFETY MEASURES IN THE PROCESS OF USE, TRANSPORT AND STORAGE

- application of safety measures: avoid excessive spreading. Good ventilation of the work area. It does not require special protective measures for the personnel who ensure the disinfection procedure.
- transport: in closed vehicles; does not require protective measures for the personnel providing the transport;
- storage: in closed original packaging, in a cool and well-ventilated place, away from heat and sunlight.

Manufacturer: **Necm Kimya**

Adress: Battalgazi Mah. Şark Cad. Aytop,
Gülöncaşı Sk. I Blok No:14/1 Sultanbeyli/İSTANBUL
Tel: +90 216 669 06 66
E-mail: info@necmkimya.com.tr

“ For Hygiene ”
BIORAD[®]



Data Sheet

READY HP WIPE

CLASSIFIED IN ACCORDANCE WITH REGULATION (EC) NO. 1272/2008 (CLP)

See the safety data sheet and label for labeling information according to Regulation (EC) No. 1272/2008 (CLP).

INFORMATION REGARDING THE ENVIRONMENT

According to the safety data sheet, the substances contained in the product comply with the biodegradability criterion established by Regulation (EC) No. 648/2004 on detergents.

Necm Kimya manufactures its products according to advanced, safe and environmentally friendly processes, economically and in compliance with high quality standards (EN ISO 9001, EN ISO 13485).

CERTIFICATION

In accordance with the European Directive on medical devices 93/42/EEC Rule 15;

- “Rule 15: All devices specifically intended to be used for disinfecting, cleaning, rinsing or, as the case may be, hydrating contact lenses are in Class IIb.
- All devices specifically intended to be used to disinfect medical devices are in Class IIa.

This rule does not apply to products that are intended to clean medical devices other than contact lenses by physical action.

Since the READY HP WIPE product is used for the disinfection of surfaces and non-invasive medical devices, the product is classified in class IIa in accordance with this norm.

Manufacturer: **Necm Kimya**

Address: Battalgazi Mah. Şark Cad. Aytop,

Gülöncaşı Sk. I Blok No:14/1 Sultanbeyli/İSTANBUL

Tel: +90 216 669 06 66

E-mail: info@necmkimya.com.tr