

MODUS FOCUSED ESWT DERMO



FOCUSED SHOCKWAVE THERAPY





Touch Screen Interface
Modus Focused ESWT - DERMO
Focused Shock Waves

MODUS FOCUSED ESWT - DERMO

ESWT (Extracorporeal Shock Wave Therapy), is considered a highly effective treatment method in the field of dermatology. This approach works by applying focused shock waves to the skin to treat skin problems. One of the fundamental principles that ESWT focuses on is stimulating skin cells and tissue structures to activate the healing mechanisms.

The effectiveness of ESWT has been proven in the treatment of serious dermatological conditions such as chronic wounds, diabetic foot ulcers, venous ulcers, and other skin lesions. The shock waves' ability to increase blood circulation and promote the growth of new blood vessels contributes to the faster healing of these types of wounds.

ESWT uses high-intensity shock waves focused on the depths of the skin to increase collagen production and aid in the revitalization of skin tissues. This provides an effective option for aesthetic concerns such as skin tightening, wrinkle reduction, and treatment of blemishes.



> MODUS FOCUSED ESWT SHOCK WAVE THERAPY DEVICE

Modus Focused ESWT Shockwave Device is a non-invasive treatment method based on the focusing of shockwaves generated outside the body onto the desired area through a fluid-filled silicone applicator. This method enhances vascularization in the affected region, activating the body's natural healing mechanisms.

- ↓ **Wide Frequency and Power Range**
Offers impulse adjustment up to 25 levels with a frequency of 4 Hz.
- ↓ **Advanced Touch Interface**
Provides a user-friendly experience with its color touch screen. Treatment parameters can be easily monitored and adjusted in real-time during sessions.
- ↓ **Patient Tracking System**
The built-in patient record and tracking menu helps manage treatment processes more effectively.
- ↓ **Preset Treatment Protocols**
The system offers visual and written guidance to ensure ease of use. Various applicator options are available for different treatment needs.
- ↓ **Automatic and Manual Control**
The device automatically stops when the predetermined number of pulses is reached, but the user can intervene manually if desired.

Modus Focused ESWT is an innovative device designed to simplify modern treatment processes. With its advanced technology optimized to meet the needs of healthcare professionals, it makes a difference at every step.

25
LEVELS OF
POWER



Modus Focused
ESWT- Dermo
Elcek & Başlık

MODUS FOCUSED ESWT HANDPIECE FEATURES

More efficient treatments with full energy transfer and ergonomic design

Advantages for the User

- > Ensures full energy transmission to anatomical regions.
- > Ergonomic design prevents hand fatigue, providing comfort during prolonged use.
- > Lightweight and well-balanced structure allows easy maneuverability.
- > Durable construction requires no maintenance, ensuring long-lasting performance.

Provides More Energy
On Deeper Tissues



F-50 Head

Penetration Depths

- 68,50 mm
- 28,50 mm
- 12,50 mm
- 84,50 mm

Modus Focused ESWT provides high treatment efficiency with a penetration depth of up to **84.50 mm**. The powerful energy reaching deep tissues helps achieve fast and effective results.



Higher Energy Output

Delivers optimal results with highly focused energy to the treatment area.



High Energy Depth

Reaches deeper tissues with effective energy, accelerating the healing process.



Fast and Effective Treatments

Saves time with shorter session durations.



Find Healing
with Shockwave
Technology

WHAT IS SHOCKWAVE THERAPY?

The method of applying therapy using shock waves from outside the body is a commonly preferred approach in various fields such as **orthopedics, physical therapy, veterinary medicine, neurology, urology, cardiology, sports medicine, and dermatology**. The advantages provided by this device include enhancing blood circulation, collagen synthesis, and oxygenation in the application area, leading to faster tissue healing. It also results in mechanically stronger tissue. Furthermore, the non-invasive nature of this system allows the desired treatment to be performed without requiring surgical procedures.

This therapy method contributes to the angiogenesis system, which involves by aiding the existing vessels, leading to the formation of new blood vessels. It also influences vasculogenesis, which is the formation of new blood vessels that accompanies the creation of new cells, through the mechanical trauma effect it creates in the tissue. By facilitating the formation of new blood vessels and initiating the repair of damage, it initiates the healing process.

Non-invasive method that provides **MUSCLE, NERVE, JOINT, BONE, TISSUE, CELL** Regeneration.

MAIN INDICATIONS

Diabetic Foot Ulcers



Diabetic foot ulcers are non-healing wounds frequently observed in diabetic patients, which can lead to serious complications if left untreated. Diabetes, due to high blood sugar levels, negatively affects nerves and blood vessels, leading to sensory loss and circulation disorders in the feet. This condition sores on wounds to form on the feet and prolongs the healing process. Diabetic foot ulcers increase the risk of infection and can result in severe outcomes such as gangrene, infection, and amputation.

Extracorporeal Shock Wave Therapy (ESWT) has emerged in recent years as a promising option for the treatment of diabetic foot ulcers. This treatment method aims to promote the healing of damaged tissues by applying high-energy shock waves in a focused manner. The fundamental mechanism of ESWT is to stimulate tissue regeneration by increasing circulation through the application of shock waves, promoting cellular restructuring, and accelerating the healing process.

Chronic Skin Ulcer



Chronic skin ulcers are non-healing and recurring open wounds that persist for a long time. These ulcers commonly develop on the legs, feet, or areas exposed to pressure. They can arise due to various reasons such as diabetic foot ulcers, venous or arterial insufficiency, pressure ulcers, and traumatic injuries. Chronic skin ulcers can negatively impact patients' quality of life and lead to serious complications due to the risk of infection. Extracorporeal Shock Wave Therapy (ESWT) is being considered as an alternative method for the treatment of chronic skin ulcers. ESWT stimulates tissue regeneration by focusing high-energy shock waves on the targeted area. This treatment can accelerate the healing process of wounds and enhance tissue regeneration, facilitating the reformation of healthy skin tissue.

Lymphedema



Lymphedema is a condition that occurs when there is a malfunction of the lymphatic system or blockage in the lymphatic channels, leading to the accumulation of fluid in the body. This accumulation can typically cause swelling and discomfort in limbs such as the arms and legs. Lymphedema can result from the surgical removal of lymph nodes, infections of the lymphatic system, or other factors that affect the lymphatic system.

Extracorporeal Shock Wave Therapy (ESWT) focuses high-energy shock waves on the targeted area, promoting tissue regeneration and enhancing circulation. In the treatment of lymphedema, the potential goal of ESWT is to stimulate lymphatic drainage, encourage the regeneration of damaged lymphatic channels, and facilitate more effective removal of accumulated fluid.

Lipedema



Lipedema is a chronic condition characterized by symmetric accumulation of fat, particularly in the lower extremities, primarily affecting women. In lipedema, fat cells in the body grow larger than normal, and fluid retention can occur. This condition typically leads to pronounced swelling and pain, especially in the hips, thighs, knees, and ankles. Lipedema is often associated with hormonal changes and may have a genetic predisposition.

Extracorporeal Shock Wave Therapy (ESWT) has been investigated as an alternative treatment method for lipedema. ESWT focuses high-energy shock waves on the targeted area, promoting tissue regeneration and enhancing circulation. In the treatment of lipedema, the potential goal of ESWT is to help reduce fluid accumulation in fat cells and regulate circulation.

Chronic Skin Ulceration



Chronic skin ulceration is characterized by the local loss of skin integrity and is often caused by underlying systemic or regional health issues such as vascular diseases, diabetes, and neurological disorders. Ulcers of this type do not heal with physical trauma, infection, and the progression of certain chronic diseases, and can lead to serious complications. ESWT, which enhances the release of endogenous angiogenic factors from cells such as endothelial cells and fibroblasts, may promote the healing process of chronic ulcers. Shock wave therapy is being considered as a potentially promising treatment method for chronic ulcers, enhancing wound healing through mechanotransduction and immunomodulatory mechanisms.

In conclusion, chronic skin ulceration is the local loss of skin integrity caused by underlying health issues and can lead to serious complications. ESWT, by increasing the release of endogenous angiogenic factors and supporting wound healing through mechanotransduction and immunomodulatory mechanisms, is being evaluated as a promising method for the treatment of chronic ulcers.

In Burn Conditions



Burns are painful injuries that occur due to thermal, electrical, chemical, or electromagnetic radiation. Smoking and exposure to open flames are common causes of burn injuries. Burns can lead to skin damage and complications, and they can be life-threatening. Burns that result in the death of skin cells can lead to serious health issues such as dehydration, electrolyte imbalance, infections, and organ failure. The prevention and treatment of burn injuries are of vital importance.

ESWT is being researched as a promising method for burn treatment. However, it is noted that the current scientific evidence is insufficient, and further research and high-quality trials are needed. More studies are required to determine the effectiveness of ESWT in the treatment of burn patients. Therefore, it is crucial to identify appropriate and effective treatment approaches for the management of burn injuries.

Venous Ulcer



Venous ulcers are chronic lesions that occur due to problems with the veins in the legs, typically found on the lower legs, ankles, or feet. These ulcers can be painful and slow to heal. Chronic venous insufficiency is one of the common causes of venous ulcers, and it results from the inadequate functioning of the valves in the leg veins, leading to the accumulation of blood in the legs. Treatment options include compression therapy, wound care, and addressing the underlying cause. ESWT has been explored as an alternative treatment option for venous ulcers, and some studies have observed that ESWT has pain-reducing and supportive effects on healing.

The study by Taheri et al. (2021) demonstrated that the use of ESWT in conjunction with Compression Bandage (KB) in the treatment of venous ulcers resulted in less pain and greater satisfaction compared to patients who received only KB treatment. ESWT significantly improved the healing process of venous ulcers, suggesting that it could be considered as a safe and effective adjunct treatment for venous ulcer therapy.

Keloid Scars




Keloid is an excessive scar tissue formation that occurs after skin trauma, extending beyond the boundaries of the initial wound and not regressing on its own. This condition often exhibits a genetic tendency, especially in individuals with darker skin tones, and while some keloids lead to an unsightly appearance, others can cause serious discomfort such as pain and itching. Treating keloids can be challenging, but innovative methods like ESWT have emerged as effective options in recent years. ESWT triggers various biological reactions in cells through the mechanical effects of shock waves and mechanotransduction. This treatment supports wound healing and tissue regeneration by providing significant improvements in volume, height, and appearance of keloids.

> MODUS FOCUSED ESWT POWERFUL AND DEEPLY EFFECTIVE SHOCKWAVE TECHNOLOGY


Modus Focused ESWT is an innovative shockwave therapy system that operates on the electrohydraulic principle. A high-voltage electrical discharge generates a powerful acoustic shockwave in a fluid medium. These waves are directed to a narrow focal point through specially designed reflectors, ensuring maximum energy transfer to the target tissue.

This technology offers a broader focal area compared to competing piezoelectric and electromagnetic ESWT systems, enabling faster and more effective treatment.



> FOCUSED SHOCKWAVES THAT TRIGGER CELL REGENERATION ↘


Focused electrohydraulic shockwaves deliver intense mechanical force and high-energy transfer to the target tissue. This biomechanical stimulation triggers cellular regeneration, accelerates blood circulation, and activates the regenerative process. With deep and precise focusing technology, it ensures rapid healing and long-lasting therapeutic effects.



> PRESSURE CHARACTERISTICS OF FOCUSED SHOCKWAVES ↘

Shock Wave Pressure Profile: The graph illustrates the pressure variation of a shock wave over time. P_+ (positive pressure) rises sharply to a peak and then gradually decreases, transitioning into negative pressure (P_-).

Pulse Width and Rise Time: The area labeled as "Pulse width (-6dB)" indicates the effective duration of the pulse, while the "Rise time" region represents the time taken for the pressure to reach its peak level. A short rise time signifies a rapid and intense impact of the shock wave.



SHOCKWAVE THERAPY IN THE TREATMENT OF CHRONIC WOUNDS

Chronic skin lesions often arise as a result of underlying conditions such as diabetes. These types of lesions can significantly impact patients' quality of life and health. However, in recent years, Extracorporeal Shock Wave Therapy (ESWT) has emerged as a promising method in the treatment of these lesions, leading to significant results. ESWT is an effective treatment approach that significantly improves and accelerates the healing processes using shock waves applied from outside the body. Research indicates that ESWT supports wound healing and provides better functionality in the restored tissue. Especially when examined in randomized controlled trials, ESWT, when used in addition to traditional wound treatment, results in approximately a 30% reduction in the wound area and an almost twofold faster healing rate. This means that the treatment process for patients can be shorter and more effective.

The effects of ESWT are not limited to just increasing blood circulation in the wound area. At the same time, the mechanical effects triggered during treatment also encourage the formation of new capillaries in the body, leading to better nourishment and healing of the tissue. This process is supported by metabolic effects, confirmed especially by the release of important proteins like eNOS and VEGF. Furthermore, it has been observed that ESWT possesses not only physical healing benefits but also antibacterial and anti-inflammatory effects. This implies that the treatment could reduce the risk of infection in the lesion area and help control inflammation.

In conclusion, the multifaceted benefits provided by ESWT in the treatment of chronic skin lesions are considered a significant advancement in the field of medicine. This treatment approach can enhance patients' quality of life and reduce the impact of issues related to skin lesions.

MODUS FOCUSED ESWT - DERMO



BEFORE

AFTER

MODUS FOCUSED ESWT - DERMO



BEFORE

AFTER

SHOCKWAVE THERAPY IN THE TREATMENT OF LIPEDEMA

Lipedema is a condition characterized by abnormal fat accumulation in the body. It is often seen in areas such as the hips, thighs, knees, and ankles, predominantly affecting women. Lipedema is a type of fat accumulation that doesn't involve the lymphatic system and is more common in individuals with a family history. This condition can develop due to hormonal factors. Lipedema leads to disproportionate fat accumulation in the body contours and typically doesn't improve with diet or exercise.

Lipedema is a condition that predominantly affects women, often demonstrating symmetrical fat accumulation in specific areas. The role of shock waves in the treatment of this condition has been investigated. Clinical studies indicate that shock waves may be effective in lipedema treatment. The applied shock waves can assist in breaking down fat cells, increasing circulation, and supporting lymphatic drainage. This can potentially lead to volume reduction in lipedema-affected areas and improvement in skin thickness.

Patients' experiences indicate that the areas where shock waves are applied feel softer, improvements in sensory functions are noted, and a noticeable reduction in the skin is perceived. Therefore, the use of shock waves in lipedema treatment is seen as a promising option for patients. However, as with any treatment method, it is important to evaluate based on individual conditions and treatment plans.

MODUS FOCUSED ESWT - DERMO



BEFORE

AFTER

MODUS FOCUSED ESWT - DERMO




BEFORE

AFTER

TECHNICAL SPECIFICATIONS

Manufacturer	İnceler Medikal Sağlık Hizmetleri San. Tic. Ltd.
Model	Modus Focused ESWT
Quality and Electrical Safety Classification	According to EN 60601-1 Class I Type B According to EN 93/42 MDD Class IIb FDA Registered Manufacturer IEC 60601-1 IEC 60601-1-2
Working Principle	Electrohydraulic
User Modes	Continuous, Burst, Auto
Treatment Start/ Stop	Main Unit Button, Foot Pedal (Optional)
Dimensions	116 mm x 387 mm x 316 mm (Main Unit) 450 mm x 350 mm x 930 mm (Trolley including)
Weight	55 kg
Voltage & Frequency	110 – 240 ±% 10 VAC & 50/60 Hz
Frequency	1 – 4 Hz
Memory	3 User Defined Treatment Protocol
Power Level	1-25
Display	Touch Screen
Electrode Lifespan	50K - 70K Shock
Operating Environment	10° C ≤ Temperature ≤ 40° C 30% Rh ≤ Humidity ≤ 80% Rh
Storage Environment	-10° C ≤ Temperature ≤ 50° C 20% Rh ≤ Humidity ≤ 90% Rh



MODUS FOCUSED ESWT
Dermo Painless and
Effective Skin Regeneration!

4 STEPS TO LIPOEDEMA TREATMENT



1. EXAMINE
Carefully inspect the area for swelling, pain, or wounds before treatment to determine suitable areas for therapy.



2. MARK
Mark the area where shock waves will be applied to ensure accurate targeting during treatment.




3. APPLY GEL
Apply medical gel to the treatment area to facilitate the effective transmission of shock waves through the skin and underlying tissues.



4. APPLY SHOCK WAVES
Firmly place the device applicator on the treatment area to deliver shock waves into the tissue. In lipedema treatment, target fat tissues to enhance circulation; in wound treatment, focus on the wound perimeter to promote tissue regeneration.



Dermatology



Technology For Health

WWW.INCELERMEDIKAL.COM



İnceler Medikal Ltd.
İvedik OSB Mah. 1472. Cad.
No: 120 Yenimahalle
Ankara/Türkiye

tel: +90 312 255 3346
sales@incelermedikal.com