

NeutrArt[®]

NEUTRAL DISPLACEMENT

Needle-Free Connector





A1010N

NeutrArt® Needle-free Valve Neutral

Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min



A1010N-PH

NeutrArt® Needle-free Valve Neutral, Photoprotective

Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min



AU1010N

NeutrArt® Needle-free Valve Neutral with Extension Line

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min



AU1010N-PH

NeutrArt® Needle-free Valve Neutral with Extension Line, Photoprotective

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min



AUWB1010N

NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line


Priming Volume	Length	Flow Rate
0.46 ml	10 cm	95 ml/min



AUWB1010N-PH

NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line, Photoprotective


Priming Volume	Length	Flow Rate
0.46 ml	10 cm	91 ml/min



A2010N

NeutrArt® Double Lumen Needle-free Valve Neutral


Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min



A2010N-PH

NeutrArt® Double Lumen Needle-free Valve Neutral, Photoprotective


Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min



A3010N

NeutrArt® Triple Lumen Needle-free Valve Neutral


Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min



A3010N-PH

NeutrArt® Triple Lumen Needle-free Valve Neutral, Photoprotective


Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min



A2010NWOR

NeutrArt® Double Lumen Needle-free Valve Neutral, Luer Slip


Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min



A2010NWOR-PH

NeutrArt® Double Lumen Needle-free Valve Neutral, Luer Slip, Photoprotective


Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min



A3010NWOR

NeutrArt® Triple Lumen Needle-free Valve Neutral, Luer Slip

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min



A3010NWOR-PH

NeutrArt® Triple Lumen Needle-free Valve Neutral, Luer Slip, Photoprotective

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min



ATU3010N

NeutrArt® Needle-free Valve with Three Way Stopcock, Neutral

Priming Volume	Length	Flow Rate
0.48 ml	15.5 cm	60 ml/min

ASSET® MEDICAL
 İkitelli OSB Mahallesi, 17. Cadde, No:17
 Basaksehir / ISTANBUL TURKEY
 T: +90 212 494 27 27 F: +90 212 494 34 37
 info@assetmedical.com - www.assetmedical.com

