

ULTRASOUND GEL

1000 ml

Ready To Use / Gel



Product Specifications

* CHEMICAL COMPOSITION

Ethanediol 1%, 2-phenoxyethanol 0.5%, Polyacrylic acid 0.67%, Glycerol 0.3% ,
2,2 -nitrioltriethanol 0.3%

* PRODUCT SPECIFICATION

In all kinds of ultrasonography and Doppler applications, it cuts the air between the skin and the probe and ensures that the ultrasonography waves come to the device screen more clearly and uninterruptedly.

* AREAS OF USE

Usage of ultrasonic gel is quite extensive. Briefly, it is used for several purposes such as ultrasound procedures, ECG procedures, cardiac stress test; ECG electrodes, epilaion procedures etc. Active ingredients do not damage the devices- the probe owing to high viscosity and vacuum technology used in production. It does not irritate the skin. It is water-soluble. It is easy to clean. It does not contain oil and oily substances. It does not contain formaldehyde and salt. It does not have toxic effect. It is odorless and suitable for all skin types.

*

* USE

It is applied sufficiently to the skin area to be treated and the procedure begins. After the procedure, the skin area is wiped with a clean paper napkin.

Product

Item.-no
Content

Ref:5743
1

Quality Features

CE

CE 2195, 93/42 EEC, CLASS 1

ULTRASOUND GEL


1000 ml

Ready To Use / Gel



Packing

Quantity Inner Packing	1
Quantity Outer Packing	20
Quantity Per Pallet	35
Product Units Per Pallet	700
Pallet Dimensions cm	120x80
LxWxH Outer Packing in mm	390x270x303
Volume in m ³ Outer Packing	0,03191 m ³
Gross Weight Outer Packing in KG	20,50 KG
Content Wight ml = gr.	1,000 gr.



Labeling

Country of Origin	Turkey
EAN Code Inner Packing	8699856751789
EAN Code Outer Packing	
Shelf Life	Two Years

Do you need more detailed information? We are happy to help!